**General instructions for  
ST270 Spear pH electrode**

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**Specifications:**

<b>Connection:</b>	BNC, Cinch
<b>Cable length:</b>	1.0m
<b>Shaft length:</b>	120mm
<b>Shaft diameter:</b>	12mm
<b>pH range:</b>	0-14pH
<b>Zero potential:</b>	7.00±0.25pH
<b>Slope:</b>	≥ 97%
<b>Impedance:</b>	≤ 100MΩ
<b>Shaft material:</b>	Glass
<b>Reference Electrolyte:</b>	3M KCl Solution
<b>Temperature range:</b>	0-100°C
<b>Temperature probe:</b>	NTC 30kΩ
<b>Reference Diaphragm:</b>	Sleeve

**OHAUS** offer different pH electrodes for different applications. Make sure you have chosen the best fit pH electrode for your samples; and start by performing a calibration of the electrode.

ST270 is a 2-in-1 pH glass electrode. A 3-in-1 electrode also has a temperature probe in it.

**Before** using the electrode in calibration or measurement, please **note**:

1. When unpacking, check the electrode for any mechanical damage to the membrane, diaphragm or body.
2. Screw off watering cap first, then the lid, rinse electrode with pure water. Wipe electrode with a tissue. Do not rub the glass membrane, or the response time might increase.
3. Eliminate any air bubbles inside the membrane. Make the electrode inside pressure in equilibrium with the atmosphere by opening the electrode refilling aperture during measurement.

### **FAQs:**

#### **1. What's the meaning of electrode condition icon?**

After calibration, the pH meter will judge the electrode slope and offset, then display the judgment on the screen:

😊 slope > 95%

Electrode condition is very good

😐 slope : 90%-95%

Electrode condition is acceptable

☹ slope : 85%-90%

Electrode condition is not good, need re-calibration or replacing.

**2. How do you properly store the pH electrode after every use?**

A: Rinse the electrode with pure water, put the electrode into the watering cap. If the solution in the cap is polluted or is dried, clean it and fill it with 3M KCl solution. Avoid placing the electrode a dry place or deionized water for a long time.

**3. What's the life of the pH electrodes?**

A: Measured samples, storage conditions and maintenance methods all could influence the life of pH electrodes; normally a pH electrode expires after 12 months. Some sample, such as strong acid or alkali, corrosive liquid will shorten the pH electrode life. The pH electrode is consumable and is generally provided with warranty.

**4. The key operation points for ST270?**

A: During use, ST270 electrode tip generally does not break. Because of the distance between tip and junction is longer than normal, the depth must exceed this distance, to ensure valid test data.

Due to the particularity of samples, ST270 electrode must be cleaned from time to time. For example, when measuring protein medium, we should use the pepsin/HCl cleaning fluid to rinse the electrode. When measuring Oily/organic solution, electrode should be rinsed by acetone or ethanol.

**Troubleshooting**

If the electrode readings are abnormal, the following procedures should be paid attention to.

1. If the slope is below 95% it needs to be re-calibrated.
2. If the measuring depth of electrode head in the sample is low the junction must be in contact with sample.
3. Soak and rinse the electrode in a cleaning solution for 2 hours after use.

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Ohaus Corporation, all rights reserved/ todos